

# Venner DruGuardian

## Improving Patient Care and Staff Safety

The unique design of Venner DruGuardian™ security device monitors and records drug cabinet access. It acts primarily to improve medication security whilst optimising staff and patient safety.

Venner DruGuardian™ can:

- Measure and improve the security of medications for ward administration<sup>1</sup>
- Reduce opportunities for drug theft and tampering<sup>1</sup>



## Features

### Alarm

Activates when drug cabinet is left open and unattended

### Locking system

To prevent tampering

### Multiple level security

3 modes of security levels Low/Medium/High

### Rechargeable battery

### Cradle

To secure Venner DruGuardian™ to the drug cabinet

### Motion/light sensor

Detects movement and activates the alarm, if the drug cabinet is left open and unattended

### Wide angle camera

Takes date and time stamped images

## Designed to protect staff

Venner DruGuardian™ was carefully designed to support staff with their clinical practice so that trolley abandonment becomes an infrequent event.<sup>1,2</sup>

**95%**  
improvement

in behaviour, when using Venner DruGuardian™ with the alarm feature.<sup>3</sup>

**20-fold**  
improvement

with the alarm delay set at 40s, drug cabinets left open and unattended for 8 min/day (alarm OFF) reduced to 0.4min/day when Venner DruGuardian™ was activated (alarm ON).

## Venner DruGuardian™ product specification

Accessories included	Cradle, Key lock
Weight	0.7kg
Dimension (including cradle)	52mm (h) x 137mm (w) x 139mm (l)
Ingress protection	IP22
Mains supply to charger	100-240 VAC, 2A charging
Backup	Backup 4.8V, 1000mAh
Power	5VDC
Battery output	Main 3.7V, 7800mAh

Product Code	Description	Qty
100000	Venner DruGuardian™	1

### References:

1. Hodges E. et al., Nursing Times, 2015: 111(40): 20-21.
2. Mariyaselvam M. 2016: Clinical Pharmacy Congress, London April 2016.
3. Mariyaselvam M. et al. 2016: Eastern Academic Health Science Network.

Venner Medical Technologies SA  
+44 (0) 1534 872555  
info@venner.com  
www.vennermedical.com

Manufactured for:  
Venner Medical (Singapore) Pte Ltd  
35 Joo Koon Circle, Singapore 629110  
www.vennermedical.com

